


# NIV PRODUCTS



## CONTENTS

### A - NIV Interfaces

Page Number

#### NASAL MASKS AND SPARE PARTS

- FlexiFit™ 407 NIV Nasal Mask

A1-A2

#### FULL FACE MASK AND SPARE PARTS

- FlexiFit™ 431 NIV Full Face Mask

A3-A4

#### ORAL MASK AND SPARE PARTS

- Oracle™ 452 NIV Oral Mask

A5-A6

**Warning:** The components in this catalogue are intended for the Non-Invasive Ventilation mask series. Mixing these parts with any other system may not achieve optimal treatment for a patient using that system. Only use NIV parts for the correct NIV mask.

*For Ventilation Products and Spare Parts please refer to the Respiratory Humidification Products Catalogue* Code 185041030  
Code 185042153 (US & Canada)

**Please note:** Some parts and/or products listed in this catalogue may not be available in all countries.

# NIV Interfaces

## *FlexiFit™ 407 NIV Nasal Mask*

Closed or Open Loop Non-Invasive Ventilation Therapy



The *FlexiFit™ 407 NIV* Nasal Mask, developed using our anthropometric database with thousands of facial contours, provides a better fit and superior comfort while maintaining the ease of use features common to the *FlexiFit™* Series.

The *FlexiFit™* Technology provides a greater range of movement at the nasal bridge, simplifying the fitting process while improving the seal and comfort. The unique *Glider™* Strap allows for freedom of movement while maintaining a seal. The new *Stretchgear™* Headgear maximizes comfort and minimizes facial pressure during head movement.

The *FlexiFit™ 407 NIV* is designed and contoured especially to meet the goal of one size to fit all, and comes completely assembled for a complete and simple Out-of-Box Solution.

Product Code	Description	Quantity
HC407 NIV	<i>FlexiFit™ 407 NIV</i> Nasal Mask includes: <ul style="list-style-type: none"><li>● Silicone Seal (one size)</li><li>● Headgear</li><li>● Glider™ Strap</li><li>● Elbows - Non-Vented &amp; Vented</li><li>● Oxygen / Pressure Port Connector</li></ul>	1 1 1 2 1

# FlexiFit™ 407 NIV Mask Spare Parts

## Elbow Kit



Product Code	Description	Quantity
400HC210	Elbows 407 NIV Non-Vented Elbow & Swivel, Vented Elbow & Vented Non-Rebreathing Valve & Instruction Manual (Australia Block)	1 of Each
400HC216	Elbows 407 NIV Non-Vented Elbow & Swivel, Vented Elbow & Instruction Manual (Australia Only)	1 of Each

## Foam Cushion




Product Code	Description	Quantity
400HC002	Foam FlexiFit™ 407	Each

## Foam Cushion & Silicone Seal Combo



Product Code	Description	Quantity
400HC501	Single Foam & Silicone Seal FlexiFit™ 407	Each

## Oxygen/Pressure Port Connector



Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack

## Forehead Pads



Product Code	Description	Quantity
900HC423	Forehead Pads - fits FlexiFit™ 407 NIV, 431 NIV	2 Pack

## Glider



Product Code	Description	Quantity
400HC211	Glider™ Strap NIV 407 NIV	Each

## Mask Headgear



Product Code	Description	Quantity
400HC301	Stretchgear™ Headgear FlexiFit™ 407	Each

**Warning:** The components in this catalogue are intended for the Non-Invasive Ventilation mask series. Mixing these parts with any other system may not achieve optimal treatment for a patient using that system. Only use NIV parts for the correct NIV mask.

# NIV Interfaces

## **FlexiFit™ 431 NIV Full Face Mask**

Closed or Open Loop Non-Invasive Ventilation Therapy



The **FlexiFit™ 431 NIV** Full Face Mask with Under Chin Design builds on the foundation of the unique features common to the **FlexiFit™** Series.

The **FlexiFit™** Technology has been incorporated into the contoured Silicone Seal to allow a large range of nasal bridge movement and simple one step auto adjustment fitting. The unique **Glider™** Strap provides freedom of movement while maintaining a seal. The crown strap and Under Chin Design provide extra stability and maintain a fit even with jaw movement. The Easy-Release Strap allows simple removal of the mask and the “Fit and Forget” Headgear ensures no on-going adjustments are required.

The **FlexiFit™ 431 NIV** Full Face is supplied fully assembled, with all sizes of Silicone Seal in the box.

Product Code	Description	Quantity
HC431 NIV	<b>FlexiFit™ 431 NIV</b> Full Face Mask includes: <ul style="list-style-type: none"><li>● Silicone Seal (Small)</li><li>● Silicone Seal (Medium)</li><li>● Silicone Seal (Large)</li><li>● Headgear with Easy Release Cord</li><li>● Glider™ Strap</li><li>● Elbows - Non-Vented &amp; Vented</li><li>● Pressure Port Covers</li></ul>	1 1 1 1 1 2 4

# FlexiFit™ 431 NIV Mask Spare Parts

## Elbow Kit



Product Code	Description	Quantity
400HC212	Elbows 431 NIV Non-Vented Elbow & Swivel, Vented Elbow & Vented Non-Rebreathing Valve & Instruction Manual (Australia Block)	1 of Each
400HC217	Elbows 431 NIV Non-Vented Elbow & Swivel, Vented Elbow & Instruction Manual (Australia Only)	1 of Each

## Silicone Seal




Product Code	Description	Quantity
400HC112	Silicone Seal FlexiFit™ 431 Small	Each
400HC104	Silicone Seal FlexiFit™ 431 Medium	Each
400HC103	Silicone Seal FlexiFit™ 431 Large	Each

## Forehead Pads



Product Code	Description	Quantity
900HC423	Forehead Pads - fits FlexiFit™ 407 NIV, 431 NIV	2 Pack

## Pressure Port Covers



Product Code	Description	Quantity
400HC213	Pressure Port Covers FlexiFit™ 431 NIV	2 Pack

## Mask Headgear




Product Code	Description	Quantity
400HC302	FlexiFit™ 431 Headgear	Each

**Warning:** The components in this catalogue are intended for the Non-Invasive Ventilation mask series. Mixing these parts with any other system may not achieve optimal treatment for a patient using that system. Only use NIV parts for the correct NIV mask.

# Oral Mask

## **ORACLE™ 452 NIV Oral Mask**

Closed or Open Loop Non-Invasive Ventilation Therapy



The Fisher & Paykel Healthcare Oracle™ was the first ever developed oral mask, a milestone in interface design.

The **Oracle™ 452 NIV** delivers pressure through the mouth rather than the nose offering the patient yet another choice of delivery. Clinically proven to function as effectively as a nasal mask, the **Oracle™ 452 NIV** is a perfect solution for mouth breathers and patients with chronic nasal obstruction who find therapy difficult with a nasal mask.

The new generation **Oracle™ 452 NIV** has exciting refinements that improve the patient experience, for ease of use and ease of seal. A new symmetrical design means the mask is always correctly fitted while the infinitely adjustable rotating dial allows for a personalized fit to provide greater comfort and simplicity.

Product Code	Description	Quantity
HC452 NIV	<b>Oracle™ 452 NIV</b> Oral Mask includes: <ul style="list-style-type: none"><li>• SoftSeal (Small and Large)</li><li>• Headgear with Easy-Release Cord</li><li>• Storage Case</li><li>• Oxygen / Pressure Port Connector</li><li>• Elbows - Non-Vented &amp; Vented</li><li>• Nasal Plugs (Small, Medium and Large)</li></ul>	1 1 1 1 2 1

# Oracle™ 452 NIV Mask Spare Parts

## Elbow Kit



Product Code	Description	Quantity
400HC214	Elbows 452 NIV Non-Vented Elbow & Swivel, Vented Elbow & Vented Non-Rebreathing Valve & Instruction Manual (Australia Block)	1 of Each
400HC218	Elbows 452 NIV Non-Vented Elbow & Swivel, Vented Elbow & Instruction Manual (Australia Only)	1 of Each

## Silicone Seal



Product Code	Description	Quantity
400HC105	Oracle™ 452 SoftSeal Small	Each
400HC106	Oracle™ 452 SoftSeal Large	Each

## Oxygen/Pressure Port Connector



Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack

## Mask Headgear



Product Code	Description	Quantity
400HC300	Oracle™ 452 Headgear Accessory plus clips	Each

## Nasal Plugs



Product Code	Description	Quantity
400HC204	Oracle™ 452 Nose Plugs S,M,L	1 of Each
400HC207	Nasal Plugs Kit 10 of Each	10 Pack

**Warning:** The components in this catalogue are intended for the Non-Invasive Ventilation mask series. Mixing these parts with any other system may not achieve optimal treatment for a patient using that system. Only use NIV parts for the correct NIV mask.

# Contacting Fisher & Paykel Healthcare

<b>International</b>	<b>Fisher &amp; Paykel Healthcare Ltd</b> P O Box 14 348, Panmure, Auckland 1134, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Sales info: info@fphcare.co.nz
<b>United States</b>	<b>Fisher &amp; Paykel Healthcare Inc</b> Tel: +1 800 446 3908 or +1 949 470 3900 Fax: +1 949 470 3933 Email: information@fphcare.com
<b>Australia</b>	<b>Fisher &amp; Paykel Healthcare Pty Ltd</b> Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 Email: info@fphcare.com.au
<b>Austria</b>	Freephone (from within Austria only) Tel: 0800 29 31 23 Fax: 0800 29 31 22
<b>Benelux</b>	<b>Fisher &amp; Paykel Healthcare SAS</b> Tel:+33 1 6446 5201 Fax:+33 1 6446 5221 Email: info@fphcare.fr
<b>China</b>	<b>Fisher &amp; Paykel Healthcare Ltd</b> Tel:+86 20 3205 3486 Fax:+86 20 3205 3495 Email: info@fphcare.cn
<b>France</b>	<b>Fisher &amp; Paykel Healthcare SAS</b> Tel:+33 1 6446 5201 Fax:+33 1 6446 5221 Email: info@fphcare.fr
<b>Germany</b>	<b>Fisher &amp; Paykel Healthcare GmbH &amp; Co KG</b> Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 Email: info@fphcare.de
<b>India</b>	<b>Fisher &amp; Paykel Healthcare Ltd</b> Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 Email:sales.india@fphcare.com
<b>Irish Republic</b>	Tel: 1800 409 011
<b>Italy</b>	Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 Email: info@fphcare.it
<b>Japan</b>	<b>Fisher &amp; Paykel Healthcare KK</b> Tel: +81 3 3661 7205 Fax:+81 3 3661 7206
<b>Netherlands</b>	Tel: +31 78 644 0924 Fax: +31 78 644 0914 Email info@fphcare.nl
<b>Northern Ireland</b>	Tel: 0800 132 189
<b>Spain</b>	Tel: +34 902 013 346 Fax +34 902 013 379 Email: sales.spain@fphcare.com
<b>Switzerland</b>	Tel: 0800 83 47 63 Fax: 0800 83 47 54
<b>United Kingdom</b>	<b>Fisher &amp; Paykel Healthcare Ltd</b> Tel: +44 1628 626 136 Fax: +44 1628 626 146 Email: info@fphcare.co.uk

w w w . f p h c a r e . c o m